



78719

## Ammonia Solution extrapure AR, 25%

Part A

CAS : 1336-21-6

Molecular Formula : NH<sub>4</sub>OH

Molecular Weight : 35.05

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Colour (APHA)	max. 10
Assay (T)	min. 25%
Density (g/ml) @ 20°C	0.91
Non Volatile Matter	max. 0.002%
Substances reducing Permanganate (O)	max. 0.0008%
Chloride (Cl)	max. 0.00005%
Sulphate (SO <sub>4</sub> )	max. 0.0002%
Iron (Fe)	max. 0.00002%
Heavy Metals (Pb)	max. 0.000005%
Phosphate (PO <sub>4</sub> )	max. 0.00005%
Carbonate (CO <sub>3</sub> )	max. 0.002%
Sulphide (S)	max. 0.00001%
Calcium (Ca)	max. 0.0001%
Copper (Cu)	max. 0.000005%
Sodium (Na)	max. 0.0005%
Potassium (K)	max. 0.0001%
Magnesium (Mg)	max. 0.00001%
Zinc (Zn)	max. 0.00002%
Pyridine & Homologues	max. 0.0002%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2672 , IMCO Class No.:8 , Packing Group:III
HSN Code	
2500 ml	28142000 (GST 18%)
500 ml	28142000 (GST 18%)

Type of Packing

2500 ml

Glass Bottle

### Available Packages

500 ml

2500 ml

25 Ltr

---

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: [info@srchem.com](mailto:info@srchem.com), website [www.srchem.com](http://www.srchem.com)

---