

Review Date: 4-Oct-2023

Section 1 - Chemical Product and Company Identification

Product Name Ammonium Ferrous Sulphate Hexahydrate extrapure AR, ACS, 99%

Product Code 61732 **CAS No** 7783-85-9

Use for Laboratory Chemicals.

Company Name Sisco Research Laboratories Pvt. Ltd.

Address 608, B Wing, Satellite Gazebo, Andheri Ghatkopar Link Road,

Andheri (E), Mumbai - 400 099, India

Section 2 - Composition/Information on Ingredients

CAS# Chemical Name: % EINECS#

7783-83-7 Ammonium Ferric Sulphate <=100 233-382-4

Section 3 - Hazards Identification

Risk advice to man

and the environment: Irritating to eyes and skin.

Section 4 - First Aid Measures

Eyes: Rinse thoroughly with plenty of water for at least 15 minutes and consult a physician.

Skin: Wash off with soap and plenty of water. Consult a physician.

Ingestion: Never give anything by mouth to an unconscious person. Rinse mouth with water. Consult a

physician.

Inhalation: If breathed in, move person into fresh air. If not breathing give artificial respiration Consult a

physician.

Notes to Physician: Consult a physician. Show this safety data sheet to the doctor in attendance

Section 5 - Fire Fighting Measures

Suitable extinguishing media: Use water spray, alcohol-resistant foam, dry chemical or

carbon dioxide.

Special protective

equipment for fire-fighters: Wear self contained breathing apparatus for fire fighting if necessary.

Section 6 - Accidental Release Measures

Personal precautions: Use personal protective equipment. Avoid dust formation. Avoid breathing dust.

Ensure adequate ventilation.

Environmental precautions: Do not let product enter drains.

Methods for cleaning up: Pick up and arrange disposal without creating dust. Keep in suitable, closed

containers for disposal.

Section 7 - Handling and Storage



www.srlchem.com

Product Code 61732

Handling: Avoid contact with skin and eyes. Avoid formation of dust and aerosols. Provide appropriate

exhaust ventilation at places where dust is formed. Normal measures for

preventive fire protection.

Storage: Store in cool place. Keep container tightly closed in a dry and well-ventilated

place.Recommended storage temperature: 2 - 8 °C.

Section 8 - Exposure Control / Personal Protection

Personal protective equipment

Respiratory protection: Where risk assessment shows air-purifying respirators are appropriate use a dust

mask type N95 (US) or type P1 (EN 143) respirator. Use respirators and components tested and approved under appropriate government standards

such as NIOSH (US) or CEN (EU).

Hand protection: The selected protective gloves have to satisfy the pecifications of EU Directive

89/686/EEC and the standard EN 374 derived from it. Handle with gloves.

Eye protection: Safety glasses

Skin and body protection: Choose body protection according to the amount and concentration of the

dangerous substance at the work

place.

Hygiene measures: Handle in accordance with good industrial hygiene and safety practice. Wash hands

before breaks and at the end of workday.

Section 9 - Physical and Chemical Properties

Physical State: Form crystalline

Molecular Formula: H4FeNO8S2 · 12H2O

Molecular Weight: 392.13

Section 10 - Stability and Reactivity

Storage stability: Stable under recommended storage conditions.

Materials to avoid: Strong oxidizing agents

Hazardous decomposition

Products: Hazardous decomposition products formed under fire conditions. - nitrogen oxides (NOx),

Sulphur oxides, Iron oxides

Section 11 - Toxicological Information



www.srlchem.com

Product Code 61732

Acute toxicity: No data available

Irritation and corrosion: No data available

Sensitisation: No data available

Chronic exposure: IARC: No component of this product present at levels greater than or equal to 0.1%

is identified as probable, possible or confirmed human carcinogen by IARC.

Genotoxicity in vitro - rat - Ascites tumor

Cytogenetic analysis

Signs and Symptoms

of Exposure: To the best of our knowledge, the chemical, physical, and toxicological properties have not

been thoroughly investigated.

Potential Health Effects

Inhalation: May be harmful if inhaled. May cause respiratory tract irritation.

Skin: May be harmful if absorbed through skin. Causes skin irritation.

Eyes: Causes eye irritation.

Ingestion: May be harmful if swallowed. **Additional Information:** RTECS: WS5900000

Section 12 - Ecological Information

Elimination information

(persistence and degradability): No data available Ecotoxicity effects: No data available

Section 13 - Disposal Considerations

Product: Observe all federal, state, and local environmental regulations. Contact a licensed professional

waste disposal

service to dispose of this material. Dissolve or mix the material with a

combustible solvent and burn in a

chemical incinerator equipped with an afterburner and scrubber.

Contaminated packaging: Dispose of as unused product.

Section 14 - Transport Information

IATA IMO RID/ADR

Shipping Name: Not Regulated for transport (Non-Haz)

Hazard Class: UN Number: Packing Group:

Section 15 - Regulatory Information



www.srlchem.com

Product Code 61732

Labelling according to EC Directives

Hazard symbols: Xi Irritant

R-phrase(s): R36/38 Irritating to eyes and skin.

S-phrase(s): S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical

advice.S36/37 Wear suitable protective clothing and gloves.

Section 16 - Other Information

Sisco Research Laboratories Pvt. Ltd. provides the information contained herein in good faith but makes no representation as to its comprehensiveness or accuracy. This document is intended only as a guide to the appropriate precautionary handling of the material by a properly trained person using this product. Individuals receiving the information must exercise their independent judgment in determining its appropriateness for a particular purpose.