

Review Date: 4-Oct-2023

Section 1 - Chemical Product and Company Identification

Product Name Barium Hydroxide Octahydrate pure, 98%

Product Code 85499 **CAS No** 12230-71-6

Use for Laboratory Chemicals.

Company Name Sisco Research Laboratories Pvt. Ltd.

Address 608, B Wing, Satellite Gazebo, Andheri Ghatkopar Link Road,

Andheri (E), Mumbai - 400 099, India

Section 2 - Composition/Information on Ingredients

CAS# Chemical Name: % EINECS# 12230-71-6 Barium hydroxide <=100 unlisted

Section 3 - Hazards Identification

EMERGENCY OVERVIEW

Danger! Causes burns by all exposure routes. Harmful if inhaled or swallowed.

Target Organs: Kidneys, central nervous system, respiratory system, gastrointestinal system, muscles,

cardiovascular system, eyes, skin.

Potential Health Effects

Eye: Causes eye burns.

Skin: Causes skin burns. May be harmful if absorbed through the skin.

Ingestion: Harmful if swallowed. Causes gastrointestinal tract burns.

Inhalation: Harmful if inhaled. Causes chemical burns to the respiratory tract.

Chronic: No information found.

Section 4 - First Aid Measures

Eyes: Immediately flush eyes with plenty of water for at least 15 minutes, occasionally lifting the upper

and lower eyelids. Get medical aid immediately.

Skin: Get medical aid immediately. Immediately flush skin with plenty of water for at least 15 minutes

while removing contaminated clothing and shoes.

Ingestion: Do not induce vomiting. Get medical aid immediately. Call a poison control center.

Inhalation: Get medical aid immediately. Remove from exposure and move to fresh air immediately. If

breathing is difficult, give oxygen. Do not use mouth-to-mouth resuscitation if victim ingested or inhaled the substance; induce artificial respiration with the aid of a pocket mask equipped with a one-way valve or other proper

respiratory medical device.

Notes to Physician: Treat symptomatically and supportively.



www.srlchem.com

Product Code 85499

Section 5 - Fire Fighting Measures

General Information: As in any fire, wear a self-contained breathing apparatus in pressure-demand,

MSHA/NIOSH (approved or equivalent), and full protective gear.

Extinguishing Media: Use water spray, dry chemical, carbon dioxide, or chemical foam.

Section 6 - Accidental Release Measures

General Information: Use proper personal protective equipment as indicated in Section 8.

Spills/Leaks: Vacuum or sweep up material and place into a suitable disposal container. Wear a self

contained breathing apparatus and appropriate personal protection. (See section). Avoid **Exposure** Controls, Personal Protection generating dusty conditions. Provide ventilation. Do not let this chemical the

environment.

Section 7 - Handling and Storage

Handling: Minimize dust generation and accumulation. Do not get in eyes, on skin, or on clothing. Do not

ingest or inhale. Use only in a chemical fume hood.

Storage: Store in a cool, dry place. Store in a tightly closed container. Corrosives area. Store under

nitrogen.

Section 8 - Exposure Control / Personal Protection

Engineering Controls: Facilities storing or utilizing this material should be equipped with an eyewash

facility and a safety shower. Use only under a chemical fume hood.

Exposure Limits:

Personal Protective Equipment

Eyes: Wear appropriate protective eyeglasses or chemical safety goggles as described by OSHA's eye and

face protection regulations in 29 CFR 1910.133 or European Standard

EN166.

Skin: Wear appropriate protective gloves to prevent skin exposure.

Clothing: Wear appropriate protective clothing to prevent skin exposure.

Respirators: A respiratory protection program that meets OSHA's 29 CFR 1910.134 and ANSI Z88.2

requirements or

European Standard EN 149 must be followed whenever workplace conditions

warrant respirator use.

Section 9 - Physical and Chemical Properties



www.srlchem.com

Product Code 85499
Physical State: Solid

Molecular Formula: Ba(OH)2.8H2O

Molecular Weight: 315.47

Section 10 - Stability and Reactivity

Chemical Stability: Absorbs carbon dioxide from the air.

Conditions to Avoid: Incompatible materials, dust generation, exposure to air.

Incompatibilities with

Other Materials: Aluminum.

Hazardous Decomposition Products: Barium oxide.

Hazardous Polymerization: Will not occur.

Section 11 - Toxicological Information

RTECS#: CAS# 12230-71-6: None listed LD50/LC50: RTECS: Not available.

Carcinogenicity: Barium hydroxide octahydrate - Not listed as a carcinogen by ACGIH, IARC, NTP, or

CA Prop 65. Barium hydroxide anhydrous - Not listed as a carcinogen by

ACGIH, IARC, NTP, or CA Prop 65.

Other:

The toxicological properties have not been fully

investigated.

Section 12 - Ecological Information

Other: Do not empty into drains. Section 13 - Disposal Considerations

Dispose of in a manner consistent with federal, state, and local regulations.

Section 14 - Transport Information

IATA IMO RID/ADR

Shipping Name: Corrosive Solids, Corrosive Solids, Corrosive Solids,

 Toxic, N.O.S.
 Toxic, N.O.S.
 Toxic, N.O.S.

 Hazard Class:
 8(6.1)
 8(6.1)
 8(6.1)

 UN Number:
 2923
 2923
 2923

Packing Group: III III III

Section 15 - Regulatory Information



www.srlchem.com

Product Code 85499

European/International Regulations

European Labeling in Accordance with EC Directives

Hazard Symbols: C

Risk Phrases: R 20/22 Harmful by inhalation and if swallowed.

R 34 Causes burns.

Safety Phrases: S 26 In case of contact with eyes, rinse immediately with plenty of water and seek

medical advice. S 36/37/39 Wear suitable protective clothing, gloves and

eye/face protection.

S 45 In case of accident or if you feel unwell, seek medical advice

immediately (show the label where possible).

WGK (Water Danger/Protection)CAS# 12230-71-6: 1, CAS# 17194-00-2

Canada CAS# 17194-00-2 is listed on Canada's DSL List

US Federal

TSCA CAS# 12230-71-6 is not on the TSCA Inventory because it is a hydrate. It is considered to be listed

if the CASnumber for the anhydrous form in on the Inventory

(40CFR720.3(u)(2)).CAS# 17194-00-2 is listed on the TSCA Inventory.

Section 16 - Other Information

Sisco Research Laboratories Pvt. Ltd. provides the information contained herein in good faith but makes no representation as to its comprehensiveness or accuracy. This document is intended only as a guide to the appropriate precautionary handling of the material by a properly trained person using this product. Individuals receiving the information must exercise their independent judgment in determining its appropriateness for a particular purpose.