



SISCO  
RESEARCH  
LABORATORIES  
PVT. LTD.

[www.srlchem.com](http://www.srlchem.com)

## Certificate of Analysis

<b>Product</b>	<b>28482 - Benzyl Chloride extrapure, 99% - [100-44-7]</b>
<b>Batch No</b>	<b>3487193</b>
<b>Analysis Date</b>	18-Nov-2023
<b>Date of Manufacture</b>	November 2023
<b>Expiry Date/Re-Test Date</b>	November 2028
<b>Molecular Formula</b>	C6H5CH2Cl
<b>Molecular Weight</b>	126.59

Test Parameters	Standards	Actual Results
Appearance (Clarity)	Clear	Clear
Appearance (Colour)	Colourless	Colourless
Appearance (Form)	Liquid	Liquid
Assay (GC)	min. 99.0%	99.61%
Density (g/ml) @ 20°C	1.099-1.100	1.099
Refractive Index (20°C)	1.538-1.539	1.538
Boiling Range	177-181°C	Passes

FRM/82/11/2/28482 Rev No. 1 Date 02/12/2021

Prepared by Sanket

In compliance with the standard specifications of SRL

Plant Site 1 : D-88/2, MIDC, Turbhe – 400 705, New Mumbai, India

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road,  
Chakala, Andheri (E), Mumbai - 400 099. India. ☎ +91-22-4268 5800

1/1

[an ISO 9001-2015, 17025:2017, CE, GMP Certified Company](#)

This COA has been generated electronically and it is valid without signature.