

Certificate of Analysis

Product

87767 - Buffer capsules pH 4 (10 capsules)

Batch No Analysis Date Date of Manufacture **3777014** 09-Jun-2024 June 2024

Expiry Date/Re-Test Date

June 2029

Test Parameters	Standards	Actual Results
pH	4.0 +/- 0.05	4.01
	(Dissolve contents of 1 capsule	
	& dilute upto 100 ml distilled water	
	to prepare solution of pH 4.0 at 20°C)	

FRM/82/11/4/87767 Rev No. 1 Date 01/11/2015

Prepared by Dnyaneshwar

In compliance with the standard specifications of SRL

Plant Site 2 : H-4, MIDC, Taloja – 410 208, Maharashtra, India