

Certificate of Analysis

Product

36018 - Cobalt (II) Chloride Hexahydrate extrapure AR,

99% - [7791-13-1]

Batch No Analysis Date Date of Manufacture Expiry Date/Re-Test Date

Molecular Formula

Molecular Weight

8027546 20-Mar-2024 March 2024 March 2029 CoCl2.6H2O

237.93

Test Parameters	Standards	Actual Results
Appearance (Colour)	Pink to red	Pink
Appearance (Form)	Crystalline compound	Crystalline compound
Solubility (Turbidity) 10% aq. solution	Clear	Clear
Solubility (Colour) 10% aq. olution	Red	Red
Assay (Complexometric)	min. 99%	99.4%
Sulphate (SO4)	max. 0.005%	Passes
Iron (Fe)	max. 0.001%	Passes
Heavy Metals (Pb)	max. 0.001%	Passes
Nitrogen Compounds (N)	max. 0.003%	Passes

FRM/82/11/4/36018 Rev No. 4 Date 17/10/2022

Prepared by Dnyaneshwar

In compliance with the standard specifications of SRL

Plant Site 2 : H-4, MIDC, Taloja – 410 208, Maharashtra, India