

Certificate of Analysis

Product

42752 - Cobalt (II) Sulphate Heptahydrate extrapure AR,

99% - [10026-24-1]

Batch No Analysis Date Date of Manufacture Expiry Date/Re-Test Date

13-Jun-2024 June 2024 June 2029

3491101

Molecular Formula Molecular Weight CoSO4.7H2O

281.10

Test Parameters	Standards	Actual Results
Appearance (Colour)	Red	Red
Appearance (Form)	Crystalline powder	Crystalline powder
Solubility (Turbidity) 10% aq. solution	Clear	Clear
Solubility (Colour) 10% aq. solution	Red	Red
Assay (Ion exchange)	min. 99%	99.3%
Chloride (CI)	max. 0.001%	Passes
Iron (Fe)	max. 0.0005%	Passes
Heavy Metals (Pb)	max. 0.001%	Passes
Copper (Cu)	max. 0.001%	Passes
Potassium (K)	max. 0.005%	Passes
Sodium (Na)	max. 0.01%	Passes

FRM/82/11/4/42752 Rev No. 2 Date 17/10/2022

Prepared by ThakurRohit

In compliance with the standard specifications of SRL

Plant Site 2 : H-4, MIDC, Taloja – 410 208, Maharashtra, India