



20259

Cyclohexane extrapure, 99%

Part A

CAS : 110-82-7

Molecular Formula : C₆H₁₂

Molecular Weight : 84.16

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Assay (GC)	min. 99%
Density (g/ml) @ 20°C	0.777-0.779
Refractive Index (20°C)	1.425-1.426
Boiling Range	80-81°C
Freezing Point	6-7°C
Acidity (HCl)	max. 0.002%
Water (KF)	max. 0.05%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1145 , IMCO Class No.:3 , Packing Group:II
HSN Code	
25 Ltr	29021100 (GST 18%)
2500 ml	29021100 (GST 18%)
500 ml	29021100 (GST 18%)
Type of Packing	
1000 ml	Glass Bottle
500 ml	Glass Bottle

Available Packages

500 ml

2500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com