



42595

Glycerol (Glycerine) Anhydrous pure, 99%

Part A

CAS : 56-81-5

Molecular Formula : C₃H₈O₃

Molecular Weight : 92.10

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Assay (NaOH titration)	min. 99%
Density (g/ml) @ 20°C	1.255 - 1.260 g/ml
Refractive Index (20°C)	1.471 - 1.475
Water (KF)	max. 1%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
25 Ltr	29054500 (GST 18%)
2500 ml	29054500 (GST 18%)
500 ml	29054500 (GST 18%)

Available Packages

500 ml

2500 ml

25 Ltr

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com