



40918

Hydrochloric Acid 0.01M Aq. Solution

Part A

CAS : 7647-01-0

Molecular Formula : HCl

Molecular Weight : 36.46

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Concentration	0.00995-0.01005M

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	24 Months
IMDG Identification	UN No.:1789 , IMCO Class No.:8 , Packing Group:II
HSN Code	
1000 ml	38229090 (GST 12%)
500 ml	38229090 (GST 12%)
Type of Packing	
1000 ml	Glass Bottle
500 ml	Glass Bottle

Available Packages

500 ml

1000 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com