



**SISCO
RESEARCH
LABORATORIES
PVT. LTD.**

www.srlchem.com

Certificate of Analysis

Product	68030 - Lead (II) Nitrate ACS, 99% - [10099-74-8]
Batch No	2524309
Analysis Date	01-Jan-2024
Date of Manufacture	January 2024
Expiry Date/Re-Test Date	January 2029
Molecular Formula	Pb(NO ₃) ₂
Molecular Weight	331.21

Test Parameters	Standards	Actual Results
Appearance (Colour)	White	White
Appearance (Form)	Solid	Solid
Assay (EDTA Titration)	min. 99.0%	99.7%
Insoluble matter	max. 0.005%	0.001%
Chloride (Cl)	max. 0.001%	Passes
Calcium (Ca)	max. 0.005%	Passes
Copper (Cu)	max. 0.002%	Passes
Iron (Fe)	max. 0.001%	Passes
Potassium (K)	max. 0.005%	Passes
Sodium (Na)	max. 0.02%	Passes

FRM/82/11/4/68030 Rev No. 2 Date 09/09/2023

Prepared by **Prajakta**

In compliance with the standard specifications of SRL

Plant Site 2 : H-4, MIDC, Talaja – 410 208, Maharashtra, India

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road,
Chakala, Andheri (E), Mumbai - 400 099. India. ☎ +91-22-4268 5800

1/1

[an ISO 9001-2015, 17025:2017, CE, GMP Certified Company](#)

This COA has been generated electronically and it is valid without signature.