



**SISCO  
RESEARCH  
LABORATORIES  
PVT. LTD.**

[www.srlchem.com](http://www.srlchem.com)

# Certificate of Analysis

**Product** 14892 - Phenol Crystalline extrapure AR, 99.5% -  
[108-95-2]  
**Batch No** 6263749  
**Analysis Date** 13-Jun-2024  
**Date of Manufacture** June 2024  
**Expiry Date/Re-Test Date** June 2029  
**Molecular Formula** C6H6O  
**Molecular Weight** 94.11

Test Parameters	Standards	Actual Results
Appearance (Colour)	White crystalline	White Crystalline
Appearance (Form)	Solid	Solid
Solubility	5% aq soln. is clear and colourless	Clear, Colourless
Assay (GC)	min. 99.5%	99.99%
Melting Point	40 - 41°C	41°C
Chloride (Cl)	max. 0.0005%	Passes
Sulphate (SO4)	max. 0.005%	Passes
Iron (Fe)	max. 0.0001%	Passes
Heavy Metals (Pb)	max. 0.0005%	Passes
Non Volatile Matter	max. 0.01%	0.0005%
Cresols (GC)	max. 0.1%	Passes
Water (KF)	max. 0.2%	0.15%
Stabilizer	0.15% H3PO2	Present

FRM/82/11/2/14892 Rev No. 2 Date 17/02/2020

Prepared by Ashish

In compliance with the standard specifications of SRL

Plant Site 1 : D-88/2, MIDC, Turbhe – 400 705, New Mumbai, India

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road,  
Chakala, Andheri (E), Mumbai - 400 099. India. ☎ +91-22-4268 5800

1/1

[an ISO 9001-2015, 17025:2017, CE, GMP Certified Company](#)

This COA has been generated electronically and it is valid without signature.