

Certificate of Analysis

Product

84749 - Potassium Hydroxide pellets extrapure, 85% -

[1310-58-3]

Batch No Analysis Date Date of Manufacture Expiry Date/Re-Test Date

21-May-2024 May 2024 May 2029

5712222

Molecular Formula Molecular Weight KOH 56.11

Test Parameters	Standards	Actual Results
Appearance (Colour)	White	White
Appearance (Form)	Hygroscopic pellets	Hygroscopic pellets
Solubility (Turbidity) 10% aq.	Clear	Clear
Solubility (Colour) 10% aq.	Colourless	Colourless
solution		
Assay (HCl titration)	min. 85%	85.3%
Chloride (CI)	max. 0.03%	Passes
Sulphate (SO4)	max. 0.02%	Passes
Carbonate (K2CO3)	max. 2%	0.2%
Iron (Fe)	max. 0.001%	Passes
Heavy Metals (Pb)	max. 0.001%	Passes

FRM/82/11/4/84749 Rev No. 2 Date 04/08/2022

Prepared by Priti

In compliance with the standard specifications of SRL

Plant Site 2 : H-4, MIDC, Taloja – 410 208, Maharashtra, India