



62673

Formic Acid extrapure AR, 98%

Part A

CAS : 64-18-6

Molecular Formula : CH₂O₂

Molecular Weight : 46.03

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Colour (APHA)	max. 10
Assay (T)	min. 98%
Density (g/ml) @ 20°C	1.210-1.220
Refractive Index (20°C)	1.370-1.371
Boiling Range	100-101°C
Freezing Point	5.5°C
Non Volatile Matter	max. 0.002%
Chloride (Cl)	max. 0.0005%
Sulphate (SO ₄)	max. 0.001%
Iron (Fe)	max. 0.0001%
Heavy Metals (Pb)	max. 0.0001%
Acetate	max. 0.3%
Sulphite (SO ₃)	max. 0.001%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1779 , IMCO Class No.:8(3) , Packing Group:II
HSN Code	
500 ml	29151100 (GST 18%)
Type of Packing	
500 ml	Glass Bottle

Available Packages

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com