



**SISCO
RESEARCH
LABORATORIES
PVT. LTD.**

www.srlchem.com

Certificate of Analysis

Product	62673 - Formic Acid extrapure AR, 98% - [64-18-6]
Batch No	1502438
Analysis Date	04-Jul-2024
Date of Manufacture	July 2024
Expiry Date/Re-Test Date	July 2029
Molecular Formula	CH2O2
Molecular Weight	46.03

Test Parameters	Standards	Actual Results
Appearance (Clarity)	Clear	Clear
Appearance (Colour)	Colourless	Colourless
Appearance (Form)	Liquid	Liquid
Colour (APHA)	max. 10	Passes
Assay (T)	min. 98%	98.74%
Density (g/ml) @ 20°C	1.210-1.220	1.210
Refractive Index (20°C)	1.370-1.371	1.370
Boiling Range	100-101°C	100-101°C
Freezing Point	5.5°C	5.5°C
Non Volatile Matter	max. 0.002%	0.0005%
Chloride (Cl)	max. 0.0005%	Passes
Sulphate (SO4)	max. 0.001%	Passes
Iron (Fe)	max. 0.0001%	Passes
Heavy Metals (Pb)	max. 0.0001%	Passes
Acetate	max. 0.3%	Passes
Sulphite (SO3)	max. 0.001%	Passes

FRM/82/11/2/62673 Rev No. 1 Date 01/11/2015

Prepared by Ashish

In compliance with the standard specifications of SRL

Plant Site 1 : D-88/2, MIDC, Turbhe – 400 705, New Mumbai, India

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road,
Chakala, Andheri (E), Mumbai - 400 099. India. ☎ +91-22-4268 5800

1/1

This COA has been generated electronically and it is valid without signature.

[an ISO 9001-2015, 17025:2017, CE, GMP Certified Company](#)