

## **Certificate of Analysis**

**Product** 

Batch No Analysis Date

**Date of Manufacture** 

Expiry Date/Re-Test Date Molecular Formula Molecular Weight 27179 - Thymol Blue - [76-61-9]

1064713

11-Feb-2024 February 2024

February 2029

C27H30O5S 466.60

Test Parameters	Standards	Actual Results
Appearance (Colour)	Brownish green	Brownish green
Appearance (Form)	Powder	Powder
Absorbance (A) of 1% aq.		
solution (pH 9.5) in a 1cm cell		
v/s H2O		
@595nm	450 - 550	542
Loss on drying	max 3%	1.96%
Suitability as indicator	Passes test	Passes

FRM/82/11/4/27179 Rev No. 1 Date 11/16/2017

Prepared by Prajakta

In compliance with the standard specifications of SRL