



CERTIFICATE OF ANALYSIS

Product Name : ZINC METAL GRANULAR 99.5% Extra Pure
Arsenic free.

Lot No. : Sample COA Analyzed On : 06-Jul-2024

Mol. Formula : Zn Mol. Weight : 65.37

Code No. : 06536 CAS No. : 7440-66-6

Mfg. Date : Jul-2024 Exp. Date : Jun-2029

HAZ. / P.G. : 9 / III UN No. : 3077

Sr.	Tests	Specifications	Results
1	Appearance	Silver gray metallic granules	Silver gray metallic granules
2	Assay	Min 99.5%	99.55%
3	Arsenic (As)	Max 0.00001%	<0.00001%
4	Iron (Fe)	Max 0.03%	0.0168%
5	Lead (Pb)	Max 0.01%	0.0014%
6	Cadmium (Cd)	Max 0.005%	0.00015%
7	Tin (Sn)	Max 0.005%	Not detectable

CONCLUSION -This above product complies as per the specifications of **LOBA CHEMIE PVT. LTD.**

Narayan
Chemist

Jagruti B. Patil
QC Manager

LOBA CHEMIE PVT. LTD.

Works : Plot No. D-22, MIDC, Tarapur Industrial Area, Tarapur, Boisar, Taluka- Palghar, Dist. Palghar, Pin-401506

Tel: 91.02525-663630/39/34

Regd Office : 107 Wode House Road, Jehanghir Villa, Colaba, Mumbai-400005

Tel: 91.22.6663 6663, Fax: 91.22.22151099

info@lobachemie.com | www.lobachemie.com